



The most reliable LC-EC applications for Drugs & Pharmaceuticals analysis

#### Antipsychotic drugs

Clozapine  
Olanzapine  
Risperidone

#### PET imaging tracer

Fluorodeoxyglucose (FDG)  
FDG impurities

#### Pharmaceuticals, API

Acetaminophen  
Artemether  
Artemisinin  
Dihydro- artemisinin  
Betadex sulfobutyl ether sodium  
Etoposide  
Epinephrine  
Heparin  
mesna BNP7787  
8-OH-DPAT  
Vincristine

#### Sulfides

Glutathione  
Aminothiols  
Disulfides

#### Aminoglycoside drugs

Amikacin  
Framycetin sulphate  
Gentamicin sulphate  
Kanamycin  
Netilmycin  
Neomycin sulfate  
Spectinomycin  
Lincomycin  
Tobramycin

## Heparin Sodium

- SweetSep™ AEX20 column
- SweetSep™ AAT amino acid trap column
- Method based on U.S. Pharmacopeia USP-NF 2023
- Limit of galactosamine in total hexosamine

### Summary

In this application note the performance of the new SweetSep™ AEX 20 anion-exchange column is evaluated for the analysis of galactosamine impurities in heparin based on the official 2023 USP Heparin Sodium monograph 2023 [1]. The SweetSep™ AEX 20 column is used in combination with the new SweetSep™ AAT amino trap column as a pre-column to retain amino acids which might be present in hydrolyzed heparin samples. Fast, reproducible, and sensitive analyses of galactosamine and glucosamine was established using the new columns with exceptionally good resolution between the two sugars. Overall, the method demonstrates excellent performance for the analysis of galactosamine impurities in heparin products.



# Heparin sodium

## Introduction

Heparin is a highly sulfated glycosaminoglycan widely used as an injectable anticoagulant. Pharmaceutical-grade heparin is derived from mucosal tissues of slaughtered meat animals such as porcine (pig) intestines or bovine (cattle) lungs.

In March 2008 a major recall of heparin was announced in the US due to reported adverse reactions (hypotension, allergic reactions) leading to fatalities in some cases [2, 3]. Upon investigation it became evident that the heparin products were adulterated with over-sulfated chondroitin sulfate, a closely related substance which resembles heparin. As a result, in 2009 the U.S. Pharmacopoeia (USP) revised the heparin monograph to address this specific adulteration issue and ensure the safety and quality of pharmaceutical-grade heparin.


In the USP-NF 2023 monograph a method based on High Performance Anion Exchange Chromatography followed by Pulsed Amperometric Detection (HPAEC-PAD) is described to analyze the organic impurities in heparin [1, 4]. Using this approach the presence of small amounts of galactosamine impurities in hydrolyzed heparin samples can be determined with high sensitivity.

## Method

The U.S. Pharmacopoeia method to determine the amount of organic impurities is based on the acid hydrolysis of heparin into glucosamine (GlcN) residues and hexuronic acid. On the contrary, over-sulfated chondroitin sulfate in adulterated heparin samples consists of galactosamine (GalN) moieties and hexuronic acid which will also be released upon hydrolysis. Both galactosamine and glucosamine can be detected by pulsed amperometric detection on an Au working electrode.



**Figure 1.** Left: SenCell with Au working electrode and Pd/H<sub>2</sub> (HyREF) reference electrode. Right: DECADE Elite electrochemical detector.



**Figure 2.** 10  $\mu$ L injection of an acid-hydrolyzed standard solution of 8  $\mu$ g/mL glucosamine and 80 ng/mL galactosamine in 50 mM HCl (hydrolyzed standard solution as described in the USP monograph).

The presence of galactosamine is a measure of the degree of contamination of heparin with over-sulfated chondroitin sulfate. The USP acceptance criteria for heparin is that not more than 1% galactosamine is present relative to the total amount of hexosamine (GlcN & GalN) in a hydrolyzed sample solution.

**Table 1**

### LC-ECD conditions

HPLC	Quaternary HPLC system
Detector	DECADE Elite electrochemical detector
Columns	SweetSep™ AEX20, 4 × 200 mm column, 5 $\mu$ m SweetSep™ AAT, 4 × 50 mm amino acid trap, 5 $\mu$ m All columns: Antec Scientific
Mobile phase	A: 14 mM KOH (elution) B: 100 mM KOH (column clean-up & regeneration) Eluents prepared & blanketed with Nitrogen 5.0
Flow rate	0.7 mL/min
System backpressure	About 240 bar
Temperature	30 °C for separation and detection
Injection volume	10 $\mu$ L
Pump piston wash	DI water (refresh weekly)
Flow cell	SenCell™ with 2 mm Au and HyREF (Pd/H <sub>2</sub> ), AST pos. 2
Potential waveform (4-step)	E1, E2, E3, E4: +0.1, -2.0, +0.6, -0.1 V ts, t1, t2, t3, t4: 0.2, 0.4, 0.02, 0.01, 0.07 s
Range	5 $\mu$ A/V (standard and sample measurement) or 500 nA/V (LOD determination)
ADF	0.5 Hz
I-cell	About 0.3 $\mu$ A



Separation

The separation of GlcN and GalN is achieved using a quaternary HPLC system equipped with a column thermostat, ET210 eluent tray for nitrogen blanketing, and DECADE Elite electrochemical detector (Figure 1). The separation is based on elution on an anion-exchange column with an alkaline mobile phase (14 mM potassium hydroxide). The analysis is based on a step-gradient, consisting of elution step, followed by a column clean-up and equilibration step (Table 2).

Table 2

Step- gradient program		
Time (min)	Mobile phase	Description
0 - 10	14 mM KOH (A)	Elution & detection
10 - 20	100 mM KOH (B)	Column clean-up/regeneration
20 - 40	14 mM KOH (A)	Equilibration to the starting condition

The eluents were carefully prepared manually using a commercial 45% KOH solution (< 0.3% K<sub>2</sub>CO<sub>3</sub>). The diluent was deionized water (resistivity >18 MΩ-cm), which was sonicated and sparged with nitrogen 5.0 prior to use. The mobile phases were blanketed with nitrogen during the analysis to minimize the build-up of carbonate ions in the mobile phase and to ensure a reproducible analysis. Take into account that the choice of a specific quaternary HPLC systems may influences the separation performance, and may require some minor tweaking of the conditions to achieve the results outlined in this application note.

In the USP Heparin monograph the use of the following column types are described: 30 x 3 mm ID amino acid trap column in series with a 30 x 3 mm ID guard column and a 150 x 3 mm ID cm analytical column. The guard and analytical column being packed with USP L69 phase. The L69 packing, as described by the USP, consists of an ethylvinylbenzene/divinylbenzene substrate agglomerated with quaternary amine functionalized 130 nm latex beads, about 6.5 μm in diameter. Antec Scientific has introduced an innovative new stationary phase, AEX20, based on the same polymeric substrate as the L69 packing, featuring highly monodisperse 5 μm particles, a substantially higher degree of crosslinking (80% instead of 55% ) in combination with quaternary amine exchange groups. Although, not specifically listed as USP L69 phase, the new stationary phase is evaluated in this application note to assess




Figure 3. 4-step PAD potential waveform for the detection of GalN and GlcN as described in the Heparin sodium USP monograph.

its performance for the separation of GalN and GlcN in hydrolyzed Heparin samples. An AEX20 analytical column with a 200 x 4 mm ID size was chosen without guard column for this evaluation. A guard column is not necessary, because of the use of an amino acid trap precolumn to retain amino acids that may be present in hydrolyzed heparin samples. The trap column will therefore also function as a guard to prevent the accumulation of contaminations and particulate matter on the analytical column.

The USP (general 621 chromatography chapter) allows adjustments of the column length (L) and particle size (dp) as long as the L/dp ratio remains in the range of –25% to + 50% of the prescribed L/dp ratio. In this case the prescribed L/dp ratio is 27700 , based on the length of the analytical + guard column (180 mm). The L/dp ratio based on the 200 mm AEX column is 40000, which is allowed as it lays within the specified upper limit of +50%. With respect to flow rate adjustments in connection with adjustments in particle size and column inner diameter the USP general chapter uses the following formula to correct the flow rate for such changes:

$$F_2 = F_1 \times [dc_2^2 \times dp_1]/(dc_1^2 \times dp_2)]$$

Where F<sub>1</sub> and F<sub>2</sub> are the flow rates for the original and modified conditions, respectively; dc<sub>1</sub> and dc<sub>2</sub> are he respective column diameters, and dp<sub>1</sub> and dp<sub>2</sub> the particle sizes. In addition, the flow rate can be adjusted by ± 50% (isocratic only).

The corrected flow rate based on the change in particle size, column length and inner diameter was calculated using the beforementioned formula. To compensate for the applied adjustments a flow rate of 0.62 mL/min was calculated.



## Heparin sodium

However, a slightly higher flow rate of 0.7 mL/min (+13%) was chosen for the AEX20 column. This flow rate is optimal for this specific column with respect to plate height (van Deemter), and well within the allowed limit for flow rate changes ( $\pm 50\%$ ) described in the USP general chromatography chapter.

### Detection

For the PAD detection of the hexoamines the Antec SenCell™ electrochemical flow cell (Figure 1) is used [5]. This flow cell with wall-jet design consists of a Au working electrode, palladium hydrogen (HyREF) reference electrode, and stainless steel auxiliary electrode. A 4-step potential waveform is used as described in the USP monograph to detect the hexoamines on the Au working electrode, see Table 1 and Figure 3. The cell current was typically about 0.3  $\mu$ A with these PAD settings. This particular 4-step waveform with a pulse duration of 500 ms has several benefits: (1) a consistent long-term peak area response and (2) minimal electrode wear [6]. The temperature for separation and detection was set to 30°C.

**Table 3**

#### USP system suitability requirements

Parameter	USP criteria	Measured
Resolution between GalN and GlcN	> 2.0	3.7
Column efficiency (GlcN)	>2000	7973
Tailing factor (GalN)	0.8-2.0	1.00
Tailing factor (GlcN)	0.8-2.0	0.98

**Table 4**

#### Linearity

Compound	Concentration range ( $\mu$ g/mL)	R
Galactosamine	0.1 - 8	0.9986
Glucosamine	0.1 - 8	0.9993

### Sample preparation

Sample digestion was achieved in the following way:

- Transfer 12 mg of heparin into a vial.
- Add 5 mL of 5 N HCl solution and vortex the solution.
- Hydrolyze the sample for 6 hours at 100 °C.
- Cool to ambient and dilute the sample 1:100 with water.

Preparation of the hydrolyzed standard solution:

- 1.6 mg/mL glucosamine stock solution: dissolve 160 mg glucosamine in 100 mL 5 N HCl.
- 16  $\mu$ g/mL galactosamine stock solution: dissolve 160 mg galactosamine in 10 mL 5 N HCl. Subsequently, add 100  $\mu$ L of the 16 mg/mL solution to 99.9 mL 5 N HCl.
- Mix equal volumes of the stock solutions (5 mL) to prepare the standard solution.
- Transfer 5 mL of the standard solution into a 7 mL screw cap vial.
- Hydrolyze the solution for 6 hours at 100 °C.
- Cool to ambient and dilute the sample 1:100 with water.

## Results

### System suitability

A chromatogram of a 10  $\mu$ L injection of an acid-hydrolyzed standard solution of 8  $\mu$ g/mL glucosamine and 80 ng/mL galactosamine in 50 mM HCl is shown in Figure 2. The retention times for galactosamine and glucosamine were 7.5 and 8.8 min, respectively. The system suitability is evaluated using the chromatogram obtained with the hydrolyzed standard solution. The results are listed in Table 3, and it is evident that all performance parameters are within the criteria of the USP system suitability requirements.

### Linearity, repeatability, and LOD


The linearity for both glucosamine and galactosamine were investigated in the concentration range of 0.1  $\mu$ g/mL – 8  $\mu$ g/mL, see Table 4. The method shows good linearity. The relative standard deviation (RSD) in peak area was determined for 10 replicate injections of the standard solution. The RSD in peak area was 0.7% and 0.3% for GalN and GlcN, respectively. The sensitivity of the method was excellent and a Limit of Detection (LOD) for galactosamine of 5.3 ng/mL was achieved.

### Sample analysis

As an example, a commercial sample was analyzed from Sigma Aldrich: Heparin sodium salt from porcine intestinal mucosa (p/n H4784, batch 051M1130V). The sample is abbreviated as sample 051M1130V from this point onwards. The sample was also spiked with the system suitability standard solution to ensure the identification of galactosamine, which may present in a small amount in the sample. The final concentration of the spike solution was 8  $\mu$ g/mL glucosamine and 80 ng/mL



galactosamine. The chromatograms obtained for the spiked sample and non-spiked hydrolyzed sample are shown in Figure 4.



**Figure 4.** 10 µL injection of non-spiked hydrolyzed sample 051M1130V (red trace) and 10 µL injection of spiked hydrolyzed sample 051M1130V (grey trace). Inset shows the zoom-in on galactosamine and glucosamine peaks.

**Table 5**

Limit of galactosamine for total hexosamine in heparin sample

Sample	USP limit %GalN	Measured %GalN
051M1160V	<1	0.02

The percentage of GalN in the hydrolyzed heparin sample is calculated compared to that of the hydrolyzed standard solution. The relative response ratio (GalNR) of GalN/GlcN in the hydrolyzed standard solution was calculated as:

(1)  $GalNR = (GalNB/GalNW) \times (GlcNW/GlcNB)$

where:

- GalNB = Peak area of GalN from hydrolyzed standard solution
- GalNW = Weight of GalN for the standard solution
- GlcNW = Weight of GlcN for the standard solution
- GlcNB = Peak area of GlcN from hydrolyzed standard solution

The percentage of galactosamine in the sample was calculated as:

(2)  $\%GalN = [(GalNU/GalNR)] / [(GalNU/GalNR) + GlcNU] \times 100$

where:

- GalNU = Peak area of GalN from hydrolyzed sample solution
- GalNR = Response ratio of GalN (1)
- GlcNU = Peak area of GlcN from hydrolyzed sample solution

The USP acceptance criteria is that not more than 1% galactosamine is present relative to the total amount of hexosamines in a hydrolyzed sample solution. The result for sample 051M1130V is listed in Table 5. The calculated %GalN is the average of a triplicate analysis of the heparin sample. The amount of galactosamine present in the sample 051M1130V is within the USP acceptance criteria for heparin.

References

1. Heparin Sodium, United States Pharmacopoeia (USP-NF) 2023, [https://doi.usp.org/USPNF/USPNF\\_M36690\\_03\\_01.html](https://doi.usp.org/USPNF/USPNF_M36690_03_01.html)
2. Information about Heparin, FDA web site: <https://www.fda.gov/drugs/postmarket-drug-safety-information-patients-and-providers/information-heparin>
3. Contaminant detected in heparin material of specified origin in the USA and in Germany; serious adverse events reported; recall measures initiated, World Health Organization, Alert No. 118, March 2008, [https://cdn.who.int/media/docs/default-source/pvg/drug-alerts/da118---alert\\_118\\_heparin.pdf?sfvrsn=bc6e1dea\\_4](https://cdn.who.int/media/docs/default-source/pvg/drug-alerts/da118---alert_118_heparin.pdf?sfvrsn=bc6e1dea_4)
4. W.R. LaCourse, “Pulsed Electrochemical Detection in High Performance Liquid Chromatography”, John Wiley & Sons, New York, 1ed, 1997.
5. H. Louw, H.J. Brouwer, N. Reinhoud, Electrochemical flowcell, (2016) US patent 9310330
6. R.D. Rocklin, A.P. Clarke, M. Weitzhandler, Anal. Chem, 70, (1998), 1496 1501

Conclusion

The DECADE Elite electrochemical detector and SenCell in combination with the new SweetSep™ AEX20 column provides a reliable solution for the sensitive analysis of low level of galactosamine in commercial heparin samples based on the official USP method. The use of the new SweetSep™ AEX20 column in combination with the SweetSep™ AAT Amino Acid Trap column resulted in an exceptionally good resolution between GalN & GlcN ensures that the system suitability requirements are easily met.



## Heparin sodium



**Figure 5.** Recommended instrument configuration for this application: the ALEXYS Carbohydrate Analyzer

The system consists of a P6.1L pump with integrated Solvent Switch Valve (SSV) capable of running step gradients, an AS6.1L autosampler, an ET 210 Eluent tray for nitrogen blanketing, and the DECADE Elite electrochemical detector. A CT 2.1 column oven with broad temperature range can be added optionally. The ALEXYS Carbohydrate Analyzer can be fully controlled by different Chromatography Data System (CDS) software: DataApex™ Clarity™ CDS (version 8.3 and up) or Thermo Scientific™ Chromeleon™ CDS (version 7.2 SR 5 and up).

### Ordering information

#### Detector only

176.0035B	DECADE Elite SCC electrochemical detector
116.4321	SenCell 2 mm Au HyREF

#### Recommended ALEXYS analyzer

180.0057W	ALEXYS Carbohydrates Analyzer - gradient (quaternary LPG)
116.4321	SenCell 2 mm Au HyREF
186.ATC00	CT2.1 Column Thermostat

#### Column

260.0020	SweetSep™ AEX20, 4 x 200 mm column, 5 µm
260.0040	SweetSep™ AAT, 4 x 50 mm amino acid trap, 5 µm

#### Software\*

195.0035	Clarity CDS single instr. incl LC, AS module
----------	--

\*) The ALEXYS Carbohydrates Analyzer can also be controlled under Thermo Fisher Scientific Chromeleon™ CDS. For the DECADE Elite electrochemical detector only also control drivers are available in Waters Empower™, Agilent OpenLab CDS and Agilent OpenLab CDS Chemstation Edition. Please contact Antec for more details.

*For research purpose only.* The information shown in this communication is solely to demonstrate the applicability of the ALEXYS system and DECADE Elite detector. The actual performance may be affected by factors beyond Antec's control and may be adjusted accordingly. Specifications mentioned in this application note are subject to change without further notice.

SweetSep, DECADE Elite, ALEXYS, SenCell, FlexCell and HyREF are trademarks of Antec Scientific. Clarity™ and DataApex™ are trademarks of DataApex Ltd. Chromeleon™ is a trademark of Thermo Fisher Scientific, Empower™ is a trademark of Waters corporation, OpenLAB™ and Chemstation are trademarks of Agilent Technologies, Inc. All other trademarks are the property of their respective owners.

#### Antec Scientific (USA)

info@AntecScientific.com  
www.AntecScientific.com  
T 888 572 0012

#### Antec Scientific (worldwide)

info@AntecScientific.com  
www.AntecScientific.com  
T +31 172 268888

